Bioinformatics file formats

Jessen V. Bredeson DOE Joint Genome Institute

Goals and outline

- Understand importance of standardized file formats
- Introduce you to commonly-used formats in bioinformatics
- Resources for manipulating or parsing them yourself

Why are (standardized) file formats important?

Data sharing and collaboration

File standards provide:
a common language for data sharing,
promote collaboration,
ensure data reusability,
reduce user errors

User1 User2 User2 compatibility Data Format/Protocols

Syntactic and semantic interoperability

"The capability to communicate, execute programs, or transfer data among various functional units in a manner that requires the user to have little or no knowledge of the unique characteristics of those units"¹


"the capability of different programs to exchange data via a common set of exchange formats, to read and write the same file formats, and to use the same protocols.... the lack of interoperability can be a consequence of a lack of attention to standardization during the design of a program"²

- 1. ISO/IEC 2382-01 Information Technology Vocabulary, Fundamental Terms
- 2. Gordon and Hernandez, The Official Guide to the SSCP Book

We have a specimen of interest...




We extract DNA...



Variants

- SNPs: Single-Nucleotide Polymorphisms
- Indels: Small (≤ 50 bp) insertions-deletions
- Structural variants: mid-to-large (> 50 bp) sized insertions, deletions, rearrangements

Variant-calling workflow



We have data, now what?

```
$ ls
Trex_genome.fasta
Trex_genome.annotation.gff3
SAMPLE_NoIndex_R1_001.fastq.gz
SAMPLE_NoIndex_R2_001.fastq.gz
$ head -3 Trex_genome.fasta
>Chr1
$ gunzip -c SAMPLE_NoIndex_R1_001.fastq.gz | head -4
@HS3:SAMPLE:1:1:1
```

Index the genome

```
$ bwa index Trex_genome.fasta # creates the following files:
$ ls Trex_genome.fasta*
Trex_genome.fasta
Trex_genome.fasta.bwt
Trex_genome.fasta.pac
Trex_genome.fasta.ann
Trex_genome.fasta.amb
Trex_genome.fasta.sa
$ samtools faidx Trex_genome.fasta # creates Trex_genome.fasta.fai
$ head -3 Trex_genome.fasta.fai
Chr1 217471166 141
                           100 101
Chr2 181034961 219646160 100 101
Chr3 153873357 402491612 100 101
$ samtools dict Trex_genome.fasta >Trex_genome.dict
$ head -3 Trex_genome.dict
@HD VN:1.0 SO:unsorted
@SQ SN:Chr1 LN:217471166 M5:56d95ce6647ea9087b857b1efa6d00dd
@SQ SN:Chr2 LN:181034961 M5:20852c561ea38c67aa67e6d655cfebf2
```

FASTA/Pearson

https://en.wikipedia.org/wiki/FASTA_format

>U31202.1 Human noggin (NOGGIN) gene, complete cds GAGCTCCGGCGGGTCAGCCGGACTGTCGGCTTCCCGGGGGCATCTGGGTCCGGCGGGGCACAGCCCTGGGC GCTGCCGAAGCCGCCGCCGCCCCCCGCGGGGGGGGGTACAGGCGGCTTCCCCCGGAGCCTGTGCAGCTCCA GAGAGAGTCAGTGGTTTCCATGGTGATGGAGCTGAAAGTGCAGGAAATTTAAAGGCTTGGACCCTGCGAG ACAGACAAACCGGTGCCAACGTGCGCGGACGCCGCCGCCGCCGCCGCCGCTGGAGTCCGCCGGGCAGAGC AGCGGCCGGNCGAAGAGCAGCGAGAGGAGGAGGGGGAGAGCGGCTCGTCCACGCGCCCTGCGCCGCCGCCG AGCCGCGGCCGCCTTCCCCAGTAGACCCGGGAGAGGAGTTGCGGCCAACTTGTGTGCCTTTCTTCCGCCC CGGTGGGAGCCGGCGCTGCGCGAAGGGCTCTCCCGGCGGCTCATGCTGCCGGCCCTGCCCCAGCC TCGGGTGAGCCGCCTCCGGAGAGACGGGGGGGGGCGCGGGGCGCGCGGGGCTCGGCGTGCTCCCCGGG GACGCGGGACGAAGCAGCCCCGGGCGCGCGCGCAGAGGCATGGAGCGCTGCCCCAGCCTAGGGGTCAC CGCCCGGCACCCAGCGACAACCTGCCCCTGGTGGACCTCATCGAACACCCCAGACCCTATCTTTGACCCCA AGGAAAAGGATCTGAACGAGACGCTGCTGCGCTGCTCGCTGCTCGGGGGCCACTACGACCCAGGCTTCATGGC AGGGCTTGGCCCAGGGCAAGAAGCAGCGCCTAAGCAAGAAGCTGCGGAGGAAGTTACAGATGTGGCTGTG GTCGCAGACATTCTGCCCCGTGCTGTACGCGTGGAACGACCTGGGCAGCCGCTTTTTGGCCGCGCTACGTG AAGGTGGGCAGCTGCTTCAGTAAGCGCTCGTGCTCCGTGCCCGAGGGCATGGTGTGCAAGCCGTCCAAGT $\mathsf{CCGTGCACCTCACGGTGCTGCGGTGGCGCTGTCAGCGGCGGGGGGGCCAGCGCTGCGGCTGGATTCCCAT$ CCAGTACCCCATCATTTCCGAGTGCAAGTGCTCGTGCTAGAACTCGGGGGCCCCCTGCCCGCACCCGGAC ACTTGATCCTCGAGCTC

>lcl|BC064885.2_cds_AAH64885.1_1 [gene=mtpn] [protein=myotrophin] [protein_id=AAH64885.1] ATGGGTGACAAGGAGTTCGTGTGGGCCATCAAGAACGGAGACCTGGATGCAGTGAAAGAATTCGTACTTG GGGGCGAGGATGTGAACCGGACGCTGGAGGAGGAGCCTATGCACTACGCTGCCGACTGCGGGCA GGATGAGGTCCTGGAGTTTCTTCTCTCGAAAGGAGCCAACATCAATGCTGCGGATAAACATGGCATCACC CCCCTACTATCTGCCTGCTACGAGGGCCCATCGCAAATGTGTCGAGTTGCTTTTATCTAAGGGAGCCGACA AGACGGTGAAGGGCCCAGACGGACTCAATGCTTTGGAATCTACAGACAACCAGGCTATCAAAGATTTGCT CCATTAA

FASTA/Pearson

FASTA Defline: Sequence ID + Description on same line, sequence string on the next

"greater than"
Start of record

Sequence ID

Required;
Any printable
non-whitespace
characters:
[!-~]

Whitespace only required if description present

>U31202.1 Human noggin (NOGGIN) gene, complete cds GASCTCCGGCGGGTCAGCCGGACTGTCGGCTTCCCGGGGGCATCTGGGTCCGGCGGGGCACAGCCCTGGGC **&**CTGCCGAAGCCGCCGCCGCCCCCCGCGGCGAGTACAGGCGGCTTCCCCCGGAGCCTGTGCAGCTCCA GAGAGAGTCAGTGGTTTCCATGGTGATGGAGCTGAAAGTGCAGGAAATTTAAAGGCTTGGACCCTGCGAG ACAGACAAACCGGTGCCAACGTGCGCGGACGCCGCCGCCGCCGCCGCCGCTGGAGTCCGCCGGGCAGAGC AGCGGCCGGNCGAAGAGCAGCGAGAGGAGGAGGGGGAGAGCGGCTCGTCCACGCGCCCTGCGCCGCCGCCG ← GCCCGGGAAGGCAGCGAGGAGCCGCCCCCGCGCCCCCGCGGTCGCCCTGGAGTAATTTCGGATGCCC AGCCGCGGCCGCCTTCCCCAGTAGACCCGGGAGAGGAGTTGCGGCCAACTTGTGTGCCTTTCTTCCGCCC CGGTGGGAGCCGGCGCTGCGCGAAGGGCTCTCCCGGCGGCTCATGCTGCCGGCCCTGCCCCAGCC TCGGGTGAGCCGCCTCCGGAGAGACGGGGGGGGGCGCGGGGCGCGCGGGGCTCGGCGTGCTCCCCGGG GACGCGGGACGAAGCAGCCCCGGGCGCGCGCGCAGAGGCATGGAGCGCTGCCCCAGCCTAGGGGTCAC CGCCCGGCACCCAGCGACAACCTGCCCCTGGTGGACCTCATCGAACACCCAGACCCTATCTTTGACCCCA AGGAAAAGGATCTGAACGAGACGCTGCTGCGCTGCTCGCTGCTCGGGGGCCACTACGACCCAGGCTTCATGGC AGGGCTTGGCCCAGGGCAAGAAGCAGCGCCTAAGCAAGAAGCTGCGGAGGAAGTTACAGATGTGGCTGTG GTCGCAGACATTCTGCCCCGTGCTGTACGCGTGGAACGACCTGGGCAGCCGCTTTTTGGCCGCGCTACGTG AAGGTGGCAGCTGCTTCAGTAAGCGCTCGTGCTCCGTGCCCGAGGGCATGGTGTGCAAGCCGTCCAAGT CCGTGCACCTCACGGTGCTGCGGTGGCGCTGTCAGCGGCGGGGGGGCCAGCGCTGCGGCTGGATTCCCAT CCAGTACCCCATCATTTCCGAGTGCAAGTGCTCGTGCTAGAACTCGGGGGCCCCCTGCCCGCACCCGGAC **ACTTGATCCTCGAGCTC**

>lcllBC064885.2_cds_AAH64885.1_1 [gene=mtpn] [protein=myotrophin] [prote

Description/Comment

CTGGGC
GCTCCA
AACCCC
Free-form text

FASTA Body/ Sequence string

Nucleotide, amino acid, IUPAC codes, alignment characters [-*]

Should be wrapped flush, but sometimes is not

FASTA files are best suffixed with ".fasta" or ".fa"; some tools require this.

CCATTAA

https://en.wikipedia.org/wiki/FASTQ_format

FASTQ Sequence Header: Sequence ID + Description on same line, sequence string on the next

Whitespace only required if description present

"At" symbol

Start of sequence portion of record

Sequence ID

Required;
Any printable
non-whitespace
characters
[!-~]

Description/Comment optional

GGATCTATGGCCATGTAGGGACCATCTGAAGGCAGATCAAAATTTCGCTGAGCAAATTTAGGGTCCGGGTTTGTT

\A<A<F--FF<-F-A7FAF-F---A<F---<FF-<F--7F----<-A7F-A----7FJ<-FF--<J<-7-FFF PSRR10178655.2 0:N:0:

ATAAAAAAAATTAATAATCTATTCTTTATTTAAAACTAATTTTTAAATTAATTGGTTTTTGTGGAATGGTAT

AAFFFJFJAJJJFFA<JF-7FJF<JJJJ--<<FJ-J<A7---<-FFJFJFJJAJ-F<F<<-F-7---7-<-<FF @SRR10178655.3 0:N:0:

CATTATATACGTCGCCACTCTTAATTTCCTTTTCCATAAGAGCGTATAATCTTGTAATACAATGTCTTCTCCAAC

AAAFFJJJJAFJJJJAFJF<JFAFJJJF-<FAJJFJ-F-F-<7-7-<FJJJJJF<JA-FF---<-7<F-F<-7-<
@SRR10178655.4 0:N:0:

AAGTTATTCTGCCTCTAATGCGATAACTGTAATCTTTAATTGTGTAATTTCTTTTTCACAATCTGAGCCACGCCA

AAAAAJF<-A<-7FJFJJJJFJFJJJJJ<FJ--7<FF-7-<--7-A<7FFJAJFFJJJAJ7FF-F7FA-7<-A-7-@SRR10178655.5 0:N:0:

AAFFFJFJJJJJJ<FJJJAFJJ<FFAJFJ--FJJJJ-FFF-<FFJFA-FJJ-AJ-<<-FFFAFJJJJJAJ-7---

FASTQ Sequence

Nucleotide, amino acid, IUPAC codes

Should *not* be wrapped flush

FASTQ files are best suffixed with ".fastq" or ".fq", some tools require this.

FASTQ Sequence Header: Paired-end or mate-pair reads

```
Read 1
                                                     Paired (or mated) reads
Type 1:
                       Read 2
                                                     may be interleaved into
    @SRR10178655.1/1
                                                 same file or separate
                                                     files. If in separate files,
    @SRR10178655.1/2
                                                        Read 1 and Read 2
                                                      sequences must be in
    -AAFFJJJAF<F-FFFFJJFFJJ<FFFJFFAJFJJJJ-F-<FJ7JJFJJF<F
                                                            same order.
Type 2:
    @SRR10178655.1 1:N:0:
    @SRR10178655.1 2:N:0:
   CATTTTTCCAAACATACCATGTCAA
                          Read 1
                      Read 2
```

FASTQ Qualities Header: Same as Sequence Header, or absent completely

"Plus" symbol_ Start of qualities portion of record

Qualities ID

Optional;
If present, typically same as Sequence ID; Must follow same rules

@SRR10178655.1 0:N:0:

GGATCTATGGCCATGTAGGGACCATCTGAAGGCAGATCAAAATTTCGCTGAGCAAATTTAGGGTCCGGGTTTGTT

AA-A<F--FF<-F-A7FAF-F---A<F---<FF-<F--7F----<-A7F-A----7FJ<-FF--<J<-7-FFFJ

TAAAAAAAAATTAATAATCTATTCTTTATTTAAAACTAATTTTTAAATTAATTGGTTTTTTGTGGAATGGTATT

AAFFFJFJAJJJFFA<JF-7FJF<JJJJJ--<<FJ-J<A7---<-FFJFJFJJAJ-F<F<<-F-7---7-<-<FFA@SRR10178655.3 0:N:0:

CATTATATACGTCGCCACTCTTAATTTCCTTTTCCATAAGAGCGTATAATCTTGTAATACAATGTCTTCTCCAAC

AAAFFJJJJAFJJJJAFJF</ri>
AAAFFJJJJAFJJJFcFAJJFJ-F-F-<7-7-<FJJJJF</pre>JA-FF---<-7<F-F</pre>
@SRR10178655.4 0:N:0:

AAGTTATTCTGCCTCTAATGCGATAACTGTAATCTTTAATTGTGTAATTTCTTTTTCACAATCTGAGCCACGCCA

AAAAAJF<-A<-7FJFJJJJFJFJJJJJJ<FJ--7<FF-7-<--7-A<7FFJAJFFJJJAJ7FF-F7FA-7<-A-7-@SRR10178655.5 0:N:0:

AAFFFJFJJJJJJJ<FFJJ

-FJJJJ-FFF-<FFJFA-FJJ-AJ-<<-FFFAFJJJJJAJ-7--

FASTQ Qualities

ASCII+offset encoded "Phred" scores.

Must be same length as sequence.

Should *not* be wrapped flush

 $Phred = -10 \bullet \log_{10}(P)$

P = fractional probability that the base call is wrong

ascii_char = chr(Phred + offset); Phred = ord(ascii_char) - offset

*					•					•					•				
<u>Dec</u>	H)	(Oct	Char	r	Dec	Нх	Oct	Html	Chr	Dec	Нх	Oct	Html (Chr	Dec	: Hx	Oct	Html Ch	<u> r</u>
0	0	000	NUL	(null)	32	20	040	 	Space	64	40	100	۵#6 4 ;	0	96	60	140	` ;	*
1	1	001	SOH	(start of heading)				!		65			A					a	a
2	2	002	STX	(start of text)	34	22	042	%#34;	rr	66	42	102	B	В	98	62	142	b	b
3				(end of text)				#		67			C					c	c
4				(end of transmission)				\$	-	68			D					d	
5	5	005	ENQ	(enquiry)				a#37;		69			E					e	
6	6	006	ACK	(acknowledge)				&					F					f	
7	7	007	BEL	(bell)	39	27	047	'	1	71			G					g	
8		010		(backspace)				a#40;		72			H					h	
9	9	011	TAB	(horizontal tab)	41	29	051))	73			I		105	69	151	i	i
10	A	012	LF	(NL line feed, new line)	42	2A	052	&# 4 2;	*				J		106	6A	152	j	j
11	В	013	VT	(vertical tab)				&#43;</td><td></td><td>75</td><td>4B</td><td>113</td><td>K</td><td>K</td><td>107</td><td>6B</td><td>153</td><td>k</td><td>k</td></tr><tr><td>12</td><td>С</td><td>014</td><td>FF</td><td>(NP form feed, new page)</td><td>44</td><td>20</td><td>054</td><td>@#44;</td><td></td><td>76</td><td>4C</td><td>114</td><td>L</td><td>L</td><td>108</td><td>6C</td><td>154</td><td>l</td><td>1</td></tr><tr><td>13</td><td>D</td><td>015</td><td>CR</td><td>(carriage return)</td><td>45</td><td>2D</td><td>055</td><td>a#45;</td><td>E 11</td><td>77</td><td>4D</td><td>115</td><td>M</td><td>M</td><td>109</td><td>6D</td><td>155</td><td>m</td><td>m</td></tr><tr><td>14</td><td>E</td><td>016</td><td>SO</td><td>(shift out)</td><td>46</td><td>2E</td><td>056</td><td>&#46;</td><td></td><td>78</td><td>4E</td><td>116</td><td>N</td><td>N</td><td>110</td><td>6E</td><td>156</td><td>n</td><td>n</td></tr><tr><td>15</td><td>F</td><td>017</td><td>SI</td><td>(shift in)</td><td>47</td><td>2F</td><td>057</td><td>/</td><td>/</td><td>79</td><td>4F</td><td>117</td><td>O</td><td>0</td><td>111</td><td>6F</td><td>157</td><td>o</td><td>0</td></tr><tr><td>16</td><td>10</td><td>020</td><td>DLE</td><td>(data link escape)</td><td>48</td><td>30</td><td>060</td><td>0</td><td>0</td><td>80</td><td>50</td><td>120</td><td>O;</td><td>P</td><td>112</td><td>70</td><td>160</td><td>p</td><td>p</td></tr><tr><td>17</td><td>11</td><td>021</td><td>DC1</td><td>(device control 1)</td><td>49</td><td>31</td><td>061</td><td>a#49;</td><td>1</td><td>81</td><td>51</td><td>121</td><td>Q</td><td>Q</td><td>113</td><td>71</td><td>161</td><td>q</td><td>q</td></tr><tr><td>18</td><td>12</td><td>022</td><td>DC2</td><td>(device control 2)</td><td>50</td><td>32</td><td>062</td><td>2</td><td>2</td><td>82</td><td>52</td><td>122</td><td>R</td><td>R</td><td>114</td><td>72</td><td>162</td><td>r</td><td>r</td></tr><tr><td>19</td><td>13</td><td>023</td><td>DC3</td><td>(device control 3)</td><td>51</td><td>33</td><td>063</td><td>3</td><td>3</td><td>83</td><td>53</td><td>123</td><td>S</td><td>S</td><td>115</td><td>73</td><td>163</td><td>s</td><td>s</td></tr><tr><td>20</td><td>14</td><td>024</td><td>DC4</td><td>(device control 4)</td><td>52</td><td>34</td><td>064</td><td>4</td><td>4</td><td>84</td><td>54</td><td>124</td><td>4;</td><td>Т</td><td>116</td><td>74</td><td>164</td><td>t</td><td>t</td></tr><tr><td>21</td><td>15</td><td>025</td><td>NAK</td><td>(negative acknowledge)</td><td>53</td><td>35</td><td>065</td><td>5</td><td>5</td><td>85</td><td>55</td><td>125</td><td>U</td><td>U</td><td>117</td><td>75</td><td>165</td><td>u</td><td>u</td></tr><tr><td>22</td><td>16</td><td>026</td><td>SYN</td><td>(synchronous idle)</td><td>54</td><td>36</td><td>066</td><td>4;</td><td>6</td><td>86</td><td>56</td><td>126</td><td>V</td><td>٧</td><td>118</td><td>76</td><td>166</td><td>v</td><td>v</td></tr><tr><td>23</td><td>17</td><td>027</td><td>ETB</td><td>(end of trans. block)</td><td>55</td><td>37</td><td>067</td><td>7;</td><td>7</td><td>87</td><td>57</td><td>127</td><td>W</td><td>W</td><td>119</td><td>77</td><td>167</td><td>w</td><td>w</td></tr><tr><td>24</td><td>18</td><td>030</td><td>CAN</td><td>(cancel)</td><td>56</td><td>38</td><td>070</td><td>8</td><td>8</td><td>88</td><td>58</td><td>130</td><td>X</td><td>Х</td><td>120</td><td>78</td><td>170</td><td>x</td><td>x</td></tr><tr><td>25</td><td>19</td><td>031</td><td>EM</td><td>(end of medium)</td><td>57</td><td>39</td><td>071</td><td>9</td><td>9</td><td>89</td><td>59</td><td>131</td><td>Y</td><td>Y</td><td>121</td><td>79</td><td>171</td><td>y</td><td>Y</td></tr><tr><td>26</td><td>1A</td><td>032</td><td>SUB</td><td>(substitute)</td><td>58</td><td>ЗА</td><td>072</td><td>:</td><td>:</td><td>90</td><td>5A</td><td>132</td><td>@#90;</td><td>Z</td><td>122</td><td>7A</td><td>172</td><td>z</td><td>z</td></tr><tr><td>27</td><td>1B</td><td>033</td><td>ESC</td><td>(escape)</td><td>59</td><td>ЗВ</td><td>073</td><td>;</td><td>;</td><td>91</td><td>5B</td><td>133</td><td>[</td><td>[</td><td>123</td><td>7B</td><td>173</td><td>{</td><td>{</td></tr><tr><td>28</td><td>10</td><td>034</td><td>FS</td><td>(file separator)</td><td>60</td><td>3С</td><td>074</td><td><</td><td><</td><td>92</td><td>5C</td><td>134</td><td>\</td><td>A.</td><td>124</td><td>7C</td><td>174</td><td> </td><td>1</td></tr><tr><td>29</td><td>1D</td><td>035</td><td>GS</td><td>(group separator)</td><td>61</td><td>ЗD</td><td>075</td><td>=</td><td>=</td><td>93</td><td>5D</td><td>135</td><td>@#93;</td><td>]</td><td>125</td><td>7D</td><td>175</td><td>}</td><td>}</td></tr><tr><td>30</td><td>1E</td><td>036</td><td>RS</td><td>(record separator)</td><td>62</td><td>ЗE</td><td>076</td><td>>;</td><td>></td><td>94</td><td>5E</td><td>136</td><td>@#94;</td><td>A</td><td>126</td><td>7E</td><td>176</td><td>~</td><td></td></tr><tr><td>31</td><td>1F</td><td>037</td><td>US</td><td>(unit separator)</td><td>63</td><td>3F</td><td>077</td><td>@#63;</td><td>2</td><td>95</td><td>5F</td><td>137</td><td>%#95;</td><td>_</td><td>127</td><td>7F</td><td>177</td><td></td><td>DEL</td></tr></tbody></table>											

P	Phred
1×10 ⁰	0
1×10 ⁻¹	10
1×10 ⁻²	20
1×10 ⁻³	30
1×10 ⁻⁴	40
1×10 ⁻⁵	50
1×10 ⁻⁶	60

Source: www.LookupTables.com

 $Phred = -10 \bullet \log_{10}(P)$


P = fractional probability that the base call is wrong

```
ascii_char = chr(Phred + offset);    Phred = ord(ascii_char) - offset
```

```
......
!"\#\%\&'()*+,-./0123456789:;<=>?@ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^-`abcdefghijklmnopqrstuvwxyz\{l\}\sim 1.00123456789:;<=>??@ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^-`abcdefghijklmnopqrstuvwxyz[\]~~.
33
                 64
                                        104
                                                     126
S - Sanger Phred+33, raw reads typically (0, 40)
          Solexa+64, raw reads typically (-5, 40)
X - Solexa
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)
J - Illumina 1.5+ Phred+64, raw reads typically (3, 41)
  with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)
  (Note: See discussion above).
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)
```

QA/QC'ing Illumina Reads

```
# Run FASTQ quality assessment tool and generate plots
$ fastqc --threads 2 --extract SRR10178655_1.fastq.gz SRR10178655_2.fastq.gz
# View the FastQC results in Safari (Mac only):
$ open -a Safari.app SRR10178655_1_fastqc/fastqc_report.html
```



Trim and align Reads

```
# Find adapter sequences in your reads and trim them off
  $ java -Xmx500m -jar ./Trimmomatic-0.39/trimmomatic-0.39.jar PE -phred33 \
      -summary SRR10178655.summary SRR10178655_1.fastq.gz SRR178655_2.fastq.gz \
     SRR10178655_1_passed.fastq.gz SRR10178655_1_failed.fastq.gz \
     SRR10178655_1_passed.fastq.gz SRR10178655_2_failed.fastq.gz MINLEN:100 \
     ILLUMINACLIP:./Trimmomatic-0.39/adapters/NexteraPE-PE.fa:2:30:10:2:keepBothReads
           Read Group tag read group ID SaMple name LiBrary name PLatform
  # Align the reads/to the genome /
   $ bwg mem -R '@RG\tID:SRR10178655\tSM:Trex\tLB:HAMMOND01\tPL:ILLUMINA' \
     Trex_genome.fasta SRR10178655_1_passed.fastq.gz SRR10178655_2_passed.fastq.gz | \
     samtools view -b - >SRR10178655.bam
  # Sort the read alignments by genome coordinate
  $ samtools sort -m 1g -o SRR10178655.srt.bam SRR10178655.bam
  # Index the BAM file for fast search (creates SRR10178655.srt.bam.bai)
   $ samtools index SRR10178655.srt.bam
Read group tag: '@RG' always
```

Read group ID: Must be unique

Sample name: Name of sample/individual/accession

Library name: Sequencing library name

Platform: Sequencing technology

http://samtools.github.io/hts-specs/SAMv1.pdf http://samtools.github.io/hts-specs/SAMtags.pdf

```
SAM: Sequence Alignment/Map format
```

BAM: Binary SAM


CRAM: Reference-Compressed SAM (also binary)

```
@HD VN:1.3 S0:coordinate
@SQ SN:ref|NC_001133| LN:230218
@SQ SN:ref|NC_001134| LN:813184
@SQ SN:ref|NC_001148| LN:948066
@SQ SN:ref|NC_001224| LN:85779
@PG ID:bwa PN:bwa VN:0.7.15-r1140 CL:bwa mem ...
@RG ID:SRR10178655 SM:Trex LB:HAMMOND01 PL:ILLUMINA
SRR10178655.85923 163 ref|NC_001133| 1 30 2S7M = 383 392
    ACATTACTC AAA))-*## NM:i:0 MD:i:7 AS:i:7 RG:Z:SRR10178655
SRR10178655.85923 83 ref|NC_001133| 383 60 9M = 1 -392
    ACCTCACAT 7JFFFFFAA NM:i:0 MD:Z:9 AS:i:9 RG:Z:SRR10178655
```

SAM Header: Meta information describing file format and data within. Header lines must start with "@" symbol (and read IDs must not). Tab separated. Reference IDs cannot be "*", "0", or "="; they have special meaning.

```
Header format version and sort order
                                                                          Read Group
                                                                        Almost required;
                                                      Program
                    VN:1.3 SO:coordinate
               @HD
                                                                        ID, sample name,
                                                 processing history
                    SN:ref|NC_001133| LN:230218
               @SQ
                                                                       and library names,
                                                  (with commands)
               @SQ
                    SN:ref|NC_001134|
                                      LN:813184
                                                                      sequencing platform
               @SQ
                    SN:ref|NC_001148|
                                      LN:948066
 Sequence
                    SN:ref|NC_001224|
               @SQ
                                      LN:85779
 Reference
               @PG
                                   VN:0.7.15-r1140 CL:bwa
                           PN:bwa
                    ID:bwa
sequence IDs
               @RG
                    ID: SRR10178655
                                    SM:Trex | B:HAMMOND01
and lengths;
                                       ref|NC_001133
               SRR10178655.85923 163
                                                                              392
listed in same
                                                               RG: Z: SRR10178655
                            AAA))-*##
 order as in
                                       refINC 0011331
               SRR10178655, 85923
                                                                              -392
                 ACCTCACAT
                                       NM:i:0
                                              MD: Z: 9
                                                       AS:i:9
                                                               RG:Z:SRR10178655
   FASTA
```

SAM Body: Describes mapping and alignment without the reference. Eleven required fields. Tab separated. Undefined values: "0" for numeric field, a "*" for non-numeric.



same rules as in FASTQ

offset = 33 required

Read Sequence and Qualities Auxiliary Info Tags Optional; RGs almost required

Reference Sequence **ID** and Position for other read in pair, if applicable

values mean other read in pair is upstream; positive is downstream

CIGAR AND Bitwise flag field details

Useful with samtools flags and samtools view -f -F

CIGAR operators

Op Meaning

M : Match

I : Insertion

D : Deletion

= : Sequence match

X : Sequence mismatch

N : Forward-skip query on reference (intron)

H : Query hard clipping

S : Query soft clipping

P : Padded reference

B : Backward-skip query on reference

Example:

Q: ATGACAGGACAGAT-GA^{GG}

R: ATG-CAGGCCAGATTGATA

3M 1I 10M 1D 2S describes same alignment as 3= 1I 4= 1X 5= 1D 2S but also reports mismatches

Bit Flags

n 2ⁿ Meaning
0 : 1 : Read is paired
1 : 2 : Read is part of proper pair
2 : 4 : Read is unmapped
3 : 8 : Other read in pair is unmapped
4 : 16 : Read is rev complemented
5 : 32 : Other read is rev complemented

6: 64: Read is R1 7: 128: Read is R2

8: 256: Alignment is a secondary hit

9 : 512 : Read fails QA/QC 10 : 1024 : Read is duplicate

11 : 2048 : Alignment is split/supplementary

To add or test for flags, use 2ⁿ values with bitwise operations:

Add flag(s) Test for flag(s) flags l = 2**0 flags & 1024 # correct flags l = 2**1 flags > 1024 # incorrect!! flags l = 2**6

QA/QC'ing Alignments

```
# Mark optical and PCR duplicate read pairs (reduce bias)
$ gatk MarkDuplicates --java-options '-Xmx1G' \
   -MAX_FILE_HANDLES 2000 \
   -I SRR10178655.srt.bam \
   -0 SRR10178655.srt.mdup.bam \
   -M SRR10178655.metrics
# Calculate QA/QC metrics for read quality etc.
$ samtools stats --ref-seq Trex_genome.fasta \
    SRR10178655.srt.mdup.bam >SRR10178655.stats
# Generate the plots
$ plot-bamstats -s Trex_genome.fasta >Trex_genome.gc
$ plot-bamstats -r Trex_genome.gc -p SRR10178655 SRR10178655.stats
# View the FastQC results in Safari (Mac only):
$ open -a Safari.app SRR10178655.html
```


QA/QC'ing Alignments

2x10⁷

1.5x10 ⁷

1x10⁷

5x10⁶



total: 54,968,582
filtered: 0 (0.0%)
non-primary: 377,364
duplicated: 0 (0.0%)
mapped: 54,878,503 (99.8%)
zero MQ: 3,226,009 (5.9%)
avg read length: 251

Bases

total: 13,797,114,082 mapped: 12,862,420,009

error rate: 1.21%

- \$ samtools view -b -f3 -F3852 SRR10178655.srt.mdup.bam > SRR10178655.srt.mdup.proper.bam
 \$ samtools index SRR10178655.srt.mdup.proper.bam
 \$ samtools tview SRR10178655.srt.mdup.proper.bam Trex_genome.fasta
 - n jessenbredeson ssh -Y bredeson@cori.nersc.gov 128×48 CCTCTACTTTCTACTGCCTCTGCATTAGCATAGGGAGAGAGGGGCGCACAGACAAGGTAGCCTTGCCGGCTAGCAATCCTCAGCGTACTCTACTTTCTGCTGCCTCTGCATTAGCATAGGGAGAGAGG

Call Variants

```
# Use local assembly of reads on the genome to calculate SNPs and Indels
gatk HaplotypeCaller \
    --minimum-mapping-quality 30 \
    --min-base-quality-score 30 \
    --read-validation-stringency SILENT \
    --reference Trex_genome.fasta \
    --input SRR10178655.srt.mdup.proper.bam \
    --output SRR10178655.vcf
```

http://samtools.github.io/hts-specs/VCFv4.3.pdf

VCF: Variant Call Format **BCF:** Binary VCF ##fileformat=VCFv4.2 ##FILTER=<ID=LowQual, Description="Locus is low quality"> ##FILTER=<ID=PASS, Description="Locus passes all filters"> ##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="The Phred-scaled prob. of the genotype"> ##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype"> ##FORMAT=<ID=AD, Number=R, Type=Integer, Description="Number of observation for each allele"> ##INFO=<ID=DP, Number=1, Type=Integer, Description="Total read depth at the locus"> ##contig=<ID=Chr1,length=217471166> ##contig=<ID=Chr2,length=181034961> #CHROM ID REF ALT FILTER POS QUAL INFO FORMAT Trex 8.826 LowQual GT:GQ:AD ./.:0:0,1Chr1 534 Α DP=1. A GT:GQ:AD 110:99:26,25 Chr1 1315 G 564.103 PASS DP=51 CTC CC 209.026 GT:GQ:AD 0|1:99:19,12 Chr1 369655 DP=31 GTT GT,GGT GT:GQ:AD 2|1:43:0,28,8 Chr1 672396 912.199 . DP=36

253.597

Chr1

GG

GGTATTTTTAG

2192815 .

GT:GQ:AD

DP=64

0/1:99:46,18

VCF Metadata Lines: For humans and computers. Required by most tools to pre-declare how to parse file body correctly. **fileformat Meta**FILTER Meta

```
Required on first line;
                                                                      explicitly defines soft
                         Tells tools how to interpret rest of file
                                                                    filters one expects to see
                                                                      in the FILTER column
##FILTER=<ID=LowQual, Description="Locus is low quality">
##FILTER=<ID=PASS, Description="Locus passes all filters">
##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="The Phred-scaled prob. of the genotype">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD, Number=R, Type=Integer, Description="Number of observation for each allele">
##INFO=<ID=DP, Number=1, Type=Integer, Description="Total read depth at the locus">
##contig=<ID=Chr1,length=217471166>
##contig=<ID=Chr2,length=181034961>
                                                           INFO
#CHROM
        POS
                                        QUAL
                                                 FILTER
Chr1
        534
                                        8.826
                                                 Low0ual
                                                          DP=1
                                                                  GT: GOFORMAT9Meta25
Chr1
        1315
                                        564.103
                                                 PASS
                                                           DP=51
                                        209. INFO. Meta
Chr1
        369655
                                                           DP = 31
                                                                  TExplicitly defines the
Chr1
        672396
                           GT, GGT
                                        Explicitly defines the = 36
                                                                  GT: Gtypes data to be 8
        21 contig Meta
Chr1
                                        types of Key=Value P=64
                                                                     observed in sample
     Optional, encouraged;
                                      data to be observed in
                                                                          column(s)
Describes reference sequences
                                           INFO column
```

observed in CHROM column

VCF Header Line: Defines columns, including the sample names. Required by most tools to parse file correctly; undefined fields set to "."

```
Locus ID
                                                            Locus-Level Meta
                            if applicable
Chromosome name
                                                                Information
                       e.g., DBsnp ID, etc.
Description="Locus is low quality" Key=Value pair info
cription="Locus Locus-level" liteabout the locus (and all
    and Position
Sequence IDs should
#be in contig Meta; ber=1, Type=IntegeQuality; Score "The samples: at the locus) f the genotype ">
#PositionsI1-basedber=1, Type=String Phred-scaled Genotype">
##FORMATA<ID=AD, Number=R, Type=Inteproblethat docus is umber of object allele">
        <1D=DP, Number=1, Type=Integern Of seally variant tal read depth at the locus">
##contia=kID=Chr1.lenath=2174
##cortig=\ID=Chr2, \length=181034961>
                                           QUAL
                                                                      FORMAT
#CHROM
                      REF ALT
                                                              INFO
                                                                                 Trex
                                                    FILTER
                                                                                   30:0,1
Chr1
         534
                                                    Low0ual
                                                                      GT: GO: AD
                                                                      T:GQ:AD
         1315
                                                    PASS
Chr1
                                                                                 0|1:99:19,12
                                                              DP=31 /GT:GQ:AD
Chr1
Chr1Reference and
                                    Locus-level
                                                            Sample-Level: AD
                                                                                 2 | Sample Field
                            GGTATTTSoft Filter(s)
                                                          Field Formatting
Chr Alternate Alleles GG
                                                                                 Contains sample
                              "PASS" = passes filters
                                                            Ordered list of
   Alleles observed in
                                                                                   genotype and
                               "." = no filters applied
                                                                                  associated info
                                                           fields present in
  reference sequence
                               anything else = failure
  and samples at the
                                                                                    at the locus
                                                               samples
          locus
```

VCF Loci: Tab-delimited columns. Alleles indexed from 0 (REF) to N (ALT) alleles. Genotypes represented with those indices

```
##fileformat=VCFv4.2
##FILTER=<ID=LowQual, Description="Locus is low quality">
##FILTER=<ID=PASS, Description="Locus passes all filters">
##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="The Phred-scaled prob. of the genotype">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
#Substitution locuser=R, Ty Complex Jocus iption="Number of observation for each allele">
##INFO=<ID=DP, Number=1, Type=Intelegription="Total read depth at the No=call or hard-
##contig=<ID=Chr1, \engtheletion and substitution!
                                                                        filtered genotype
##contig=<ID=Chr2,length=181034961>
                 ID
                    REF
                                       QUAL
                                                                 FORMAT
                                                                           Trex
#CHROM
        POS
                          ALT
                                                FILTER
                                                          INFO
                                       8.826
                                                LowQual
                                                                 GT:GQ:AD
                                                                           ./.:0:0,1
Chr1
        534
                          Α
                                                         DP=1
                          G
                                                                           110:99:26,25
                                                         DP=51
                                                                 GT:GQ:AD
Chr1
                                       564.103
                                                PASS
        1315
                                                                           01:99:19,12
                     CTC
                          CC
                                       209.026
                                                                 GT:GQ:AD
Chr1
        369655
                                                         DP=31
                                                                           211:43:0,28,8
                          GT, GGT
                                                                 GT:GQ:AD
        672396
Chr1
                     GTT
                                       912.199 .
                                                         DP=36
                     GG
                                                                 GT:GQ:AD 0/1:99:46,18
Chr1
        2192815
                          GGTATTTTTAG
                                       253.597
                                                          DP=64
                                               Phased genotypes
Deletion locus
                                                                             Allele Depth
```

Insertion locus

Unphased genotype

Read count for

each allele

Annotation files

BED


https://genome.ucsc.edu/FAQ/FAQformat.html#format1

BED: Browser Extensible Data format

```
track itemRgb="0n"
Chr1
     0
             126500
                    CLEAR1
                           0 + 0
                                        126500 0,0,0
                                        128500 213,221,213
     126500 128500
                    BREAK1
                           0 + 126500
Chr1
     128500
            278000 CLEAR2 0 + 128500
Chr1
                                        278000
                                               0,0,0
     278000
            280000
                    BREAK2 0 + 278000
                                               213,221,213
Chr1
                                        280000
Chr1
     280000
             362500 CLEAR3 0 + 280000
                                        362500
                                               0,0,0
     362500
             366000
                          0 + 362500
                                        366000
                                               213,221,213
Chr1
                    BREAK3
                                               0,0,0
Chr1
     366000
            427500 CLEAR4
                           0 + 366000
                                        427500
                           0 + 427500
Chr1
     427500
            429500
                    BREAK4
                                        429500
                                               213,221,213
Chr1
     429500
             599500
                    CLEAR5
                           0 + 429500
                                        599500
                                               0,0,0
Chr1
     599500
            605500
                    BREAK5
                                599500
                                        605500
                                               213,221,213
```

BED

BED: Columns tab-delimited. First three required, all others optional (first 6 typical).



GFF3

https://github.com/The-Sequence-Ontology/Specifications/blob/master/gff3.md

GFF: Generic Feature Format

```
##gff-version 3
##species http://www.ncbi.nlm.nih.gov/Taxonomy/Browser/wwwtax.cgi?id=436495
##genome-build RexBase Trex1
##sequence-region Chr1 1 217471166
# Note Trex_genome.fasta, complete genome
Chr1 Gnomon gene
                          43895 78350 . + .
                                                ID=gene32251; Name=L0C101732307
                                                ID=rna61088; Name=XM_012954515.1; Parent=gene32251
             mRNA
                          43895 78350
Chr1 Gnomon
                          43895 43947 .
Chr1 Gnomon
            CDS
                                             0 ID=rna61088.1.CDS;Parent=rna61088
Chr1 Gnomon exon
                          43895 43947
                                                ID=rna61088.1.exon;Parent=rna61088
                         43895 43897
                                             0 ID=rna61088.1.start_codon;Parent=rna61088
Chr1 Gnomon start_codon
             CDS
                          48839 49007
                                            1 ID=rna61088.2.CDS; Parent=rna61088
Chr1 Gnomon
                                          +
                          48839 49007
Chr1 Gnomon
                                             . ID=rna61088.2.exon; Parent=rna61088
             exon
                                          +
Chr1 Gnomon
             CDS
                          53889 54000
                                          + 0 ID=rna61088.3.CDS; Parent=rna61088
                                            . ID=rna61088.3.exon;Parent=rna61088
Chr1 Gnomon
                          53889 54000
             exon
                                            2 ID=rna61088.4.CDS; Parent=rna61088
                          55055 55173
             CDS
Chr1
     Gnomon
                                          +
                          55055 55173
                                                ID=rna61088.4.exon; Parent=rna61088
Chr1
     Gnomon
             exon
```

GFF3

GFF Header: Pragma begin with "##", comments with "#". Format pragma required for GFF3, highly-recommended for GFF2/GTF.

Pragma/Directives

Pre-declared set of pragma with specific formats/definitions.

Mostly for computers/browsers.

```
##qff-version 3
##spectes http://www.ncbi.nlm.nih.gov/Taxonomy/Browser/wwwtax.cgi?id=436495
##genome-build RexBase Trex1
##sequence-region Chr1 1 217471166
# Note Trex_genome.fasta, complete genome
Chr1 Gnomon gene
                           43895
                                                  ID=gene32251; Name=L0C101732307
ChrFormatiVersion
                                                  ID=rna61088; Name=XM_012954515.1; Parent=gene32251
                           43895
                           43895
                                                  ID=rna61088.1.CDS; Parent=rna61088
 Pragma/Directive
                                  Comments
                           43895
                                                  ID=rna61088.1.exon; Parent=rna61088
Required for GFF3, don
                                                  ID=rna61088.1.start_codon;Parent=rna61088
                           43895
                                 Free-form text
                                                  ID=rna61088.2.CDS; Parent=rna61088
                           48839
highly-recommended
                                  for humans,
                                                  ID=rna61088.2.exon; Parent=rna61088
                           48839
for GFF2/GTF formats
                           53889
                                  5ignored by 0
                                                  ID=rna61088.3.CDS; Parent=rna61088
                           53889
                                                  ID=rna61088.3.exon; Parent=rna61088
      Gnomon
                           55055
                                                  ID=rna61088.4.CDS; Parent=rna61088
      Gnomon
              CDS
                                                  ID=rna61088.4.exon; Parent=rna61088
      Gnomon exon
                           55055
```

GFF3

GFF Features: Nine tab-delimited fields required. Null values a "."

```
Reference ID

Chromosome/scaffold ID

May only contain

characters in set:

[a#ZA-Z0-9:^^*$@!+_?-]

##species nttp://www.ncbi.nlm.

##genome-build RexBase frex1

##sealence-region Ch 1 1 21747

# Note Trex_genome fasta, compared the compared the compared to the compared
```

Chr1

Chr1

Gnomon

Gnomon

Gnomor

Gnomon

Feature Type

Must be SO term or accession number

Score

nlm.nih.gov/Taxonomy/Bro ex1 floating point 21747116@number complete genome 43895 78350 . +

43895 78350 . + . 43895 43947 . + 0

43895 43947 . + 43895 43897 . + 6

48839 53889

Chriusually the program of 889 Chriusually the program of 55055 Chriusually the program of 55055 Chriusually the program of 55055

exon

CDS

Gnomon Source

start_codon

generated the annotations

Feature Strand

either "+" or "-", or "." if no strand applies **Feature Attributes**

Semi-colon separated Key=Value pairs; reserved keys begin with capitals letters;

"Parent" attribute defines feature hierarchy; must use

www.cgi?id=436495 URL-escaping for forbidden characters

ID=gene32251; Name=L0C101732307

ID=rna61088; Name=XM_012954515.1; Parent=gene32251

ID=rna61088.1.CDS;Parent=rna61088
ID=rna61088.1.exon;Parent=rna61088

ID=rna61088.1.start_codon;Parent=rna61088

D=rna61088.2.CDS;Parent=rna61088

D=rna61088.2.exon;Parent=rna61088 D=rna61088.3.CDS;Parent=rna61088

D=rna61088.3.exon;Parent=rna61088

ID-Codon Phase arent=rna61088

either 0, 1, expra?Parent=rna61088

Offset to next codon position

Start and End

55**Positions**² 5517³based

coordinates on "+" strand

Resources

File manipulation/filtering

pysam (API)	FASTA/Q, BED, B/CR/SAM, B/VCF	https://pysam.readthedocs.io/en/latest/api.html#sam-bam-cram-files
BioPython	Many	https://biopython.org
pyFaidx (API)	FASTA	https://doi.org/10.7287/peerj.preprints.970v1
Seqtk	FASTA/Q	https://github.com/lh3/seqtk
Seqkit	FASTA/Q	https://doi.org/10.1371/journal.pone.0163962
seqmagick	Many	https://seqmagick.readthedocs.io
bedtools	BAM, BED, GFF, VCF	https://bedtools.readthedocs.io
bcftools	B/VCF	https://samtools.github.io/bcftools
genometools	FASTA/Q, GFF, GTF	http://genometools.org
gffread & gffcompare	GFF, GTF	https://github.com/gpertea/gffread https://github.com/gpertea/gffcompare
samtools	FASTA/Q, B/SAM	https://github.com/samtools/samtools
bamtools	B/SAM	https://github.com/pezmaster31/bamtools
vcftools	B/VCF	https://vcftools.github.io/man_latest.html
Picard	FASTA/Q, BED, B/CR/SAM, B/VCF	https://broadinstitute.github.io/picard/

Resources

QA/QC, Adapter and Quality trimming

trimmomatic	FASTQ	http://usadellab.org/cms/?page=trimmomatic
FastQC	FASTQ, B/SAM	https://www.bioinformatics.babraham.ac.uk/projects/fastqc/
Sickle	FASTA/Q	https://github.com/ucdavis-bioinformatics/sickle
Scythe	FASTA/Q	https://github.com/ucdavis-bioinformatics/scythe
Sabre	FASTA/Q	https://github.com/najoshi/sabre
cutadapt	FASTA/Q	https://cutadapt.readthedocs.io/en/stable/

Alignment

minimap2	FASTA/Q	https://github.com/lh3/minimap2
miniprot	FASTA	https://github.com/lh3/miniprot
BWA	FASTA/Q	https://github.com/lh3/bwa
hisat2	FASTA/Q	https://daehwankimlab.github.io/hisat2/
STAR	FASTQ	https://github.com/alexdobin/STAR
GMAP	FASTA/Q	http://research-pub.gene.com/gmap/
exonerate	FASTA	https://www.ebi.ac.uk/about/vertebrate-genomics/software/exonerate

Resources

Variant calling

FreeBayes	BAM, VCF	https://github.com/ekg/freebayes
GATK4	FASTA/Q, B/CRAM, VCF	https://software.broadinstitute.org/gatk/documentation
DeepVariant	FASTA/Q	https://github.com/google/deepvariant
vg	FASTA/Q	https://github.com/vgteam/vg

Common file issues

- Non-printable characters
- Non-ASCII encoded characters
- Incorrect formatting (spaces instead of tabs)
- Truncated files

OC -C